

CLEANING YOUR PUFFER DEVICE




Your reliever puffer (Asmol/Ventolin/Zempreon) should be cleaned **once a week**. To do so, follow the steps below.



1

Remove the metal canister and dose counter from the plastic casing (actuator) of the puffer and remove the mouthpiece cover.



2

Rinse the actuator thoroughly under warm running water.
Please do not put the metal canister in water.



3

Dry the actuator thoroughly inside and out.



4

Replace the metal canister and the mouthpiece cover.

NOTE: If your inhaler fails to provide relief within approximately 4 minutes of use, it may have become clogged. Spray the inhaler into the air and observe the flow from the mouthpiece, clean if required.

If you need to use your inhaler before it is completely dry, shake off as much excess water as possible before using. Spray one spray out to clear out any remaining water droplets. After use, repeat cleaning steps and allow to dry.



NOT BLOCKED



BLOCKED

The manufacturer of VENTOLIN and ZEMPREON (GSK) would like people who have had any problems with the new devices with dose counter to contact GSK Medical Information on 1800 033 109. The manufacturer of Asmol is Mylan Australian please email medinfo_anz@mylan.com or call 1800 274 276. You can report concerns to the Therapeutic Goods Administration, www.tga.gov.au/reporting-problems. Feel free to contact us on **1800 ASTHMA** (1800 278 462)